

# Bougies AG50++

– Dilation of meatus before endoscopy



**AG50++ Filiform bougie, straight olive tip****Manufacturer**

COLOPLAST A/S

**Distributor**

COLOPLAST A/S

**CE class**

Is (Rule 5)

**Sterilisation**Sterile and single use.  
Ethylene oxide (EO)**Indications**Exploration of the meatus and urethra.  
*Instructions for use SH2001***Contraindication**

Any known allergies to the medical device materials.

Marketing

**References**

CH / FR	Order references
3	AG5003
4	AG5004
5	AG5005
6	AG5006
7	AG5007

**Description**Length (cm) : 30  
Distal tip : straight, olive tip  
Proximal tip : straight cut**Materials** : NEOPLEX<sup>®</sup>**Containing Latex** : No**Shelf life**

5 years

**Packaging**

<b>Use unit</b>	1 bougie in a peel pouch
<b>Order unit</b>	1 use unit in a carton box (mm) 100x555

Regulatory Affairs

Frsbu 16/06/10

