

# Bougies AG53++

– Dilation of meatus before endoscopy



**AG53++ Straight olive tip bougie****Manufacturer**

COLOPLAST A/S

**Distributor**

COLOPLAST A/S

**CE class**

Is (Rule 5)

**Sterilisation**Sterile and single use.  
Ethylene oxide (EO)**Indications**Dilation of urethral stenoses.  
*Instructions for use SH2001***Contraindication**

Any known allergies to the medical device materials.

**References**

CH / FR	Order references
6	AG5306
7	AG5307
8	AG5308
9	AG5309
10	AG5310
12	AG5312
14	AG5314
16	AG5316
18	AG5318
20	AG5320
22	AG5322
24	AG5324
26	AG5326
28	AG5328
30	AG5330

**Description**Length (cm) : 34  
Distal tip : straight olive tip  
Proximal end : straight**Materials:** NEOPLEX®**Containing Latex:** No**Shelf life**

5 years

**Packaging****Use unit** 1 peel pouch  
**Order unit** 1 use unit in a carton box (mm)  
100x555

Marketing

Regulatory Affairs



Frsbu 16/06/10

