


# Ballonkateter silikone

## Kateter med 2 øjne over ballon



## AJ66++ Short-term nephrostomy catheter with balloon in silicone

### Manufacturer

COLOPLAST A/S

### Distributor

COLOPLAST A/S

### CE class

IIb (Rule 8)

### Sterilisation

Sterile and single use  
Ethylene oxide (EO)

### Indications

Short-term drainage of the upper urinary tract by percutaneous way.

*Instructions for use SH2105*

### Shelf life

5 years

Marketing



### References

CH / FR	Length (cm)	Balloon volume (ml)	Valve colour	Code
10	30	3	Violet	AJ6610
12	43	5	White	AJ6612
14	43	5	Green	AJ6614
16	43	5	Orange	AJ6616
18	43	5	Red	AJ6618
20	43	5	Yellow	AJ6620

### Description

Distal tip: open, straight and tapered  
Proximal end: Fitted with a funnel and a valve of which the colour varies according to the size of the catheter  
Channels: 2  
Eye: 1 through and through  
Graduation

**Material:** Silicone

**Containing latex:** No

### Packaging

**Use unit** 1 catheter in a pouch and peel pouch

**Order unit** 10 uses units in a carton box (mm) 90x65x550

Regulatory Affairs

