

Kateter Latex – couvellaire

**AM36++ COUVELAIRE prostatectomy catheter
in reinforced latex****Manufacturer**

COLOPLAST A/S

Distributor

COLOPLAST A/S

CE class

IIa (Rule 5)

SterilisationSterile and single use
Ethylene oxide (EO)**Indications**

Short-term drainage of the vesical urines, after prostate surgery.
Post-operative irrigation of the bladder.
Post-operative haemostasis of the prostate area.
The data in the literature recommend that latex urinary catheters should not be left in place for more than 7 days.

Contraindications

The same as for urethral urinary catheterisation and generally, known allergic reactions due to the device material (e.g. latex) or sterilisation agents (particularly ethylene oxide and its derivatives).
CAUTION : These products contain natural latex, which may cause allergic reactions.

Instructions for use 2SH2115

Marketing

**References**

CH / FR	Ring colour	Order references
18	Red	AM3618
20	Yellow	AM3620
22	Violet	AM3622
24	Blue	AM3624

Description

Length (cm) : 42
Distal tip : straight bevel tip
Proximal end : funnel, with reinforcement coil, rubber ring and Luer valve
Channels : 3 (drainage, irrigation, inflation)
Eyes : 2 opposite and staggered
Balloon (ml) : 30 - 50
Length (cm) : 42

Materials : Reinforced latex**Shelf life**

5 years

Packaging

Use unit	1 pouch and peel pouch
Order unit	5 use units in a carton box (mm) 90x40x550

Regulatory Affairs

