

BLADDER EVACUATOR

References

Order reference	Description
000PABVAC	Bladder evacuator

Description

Device description:

Bladder evacuator, bottle with a pre-attached cap and filter and silicone tube.

Nozzle Adaptor for Storz/Wolf/Olympus/ACMIresectoscope and a caplug (Endcap).

Materials:

Bottle: EVA (Ethylene Vinyl Acetate)

Connector: EVA (Ethylene Vinyl Acetate)

Flexible tube: Silicone

Cap plug : LDPE

Bottle cover with perforated tube : PP (Polypropylene)

O-ring : Silicone

Not manufactured with Phthalates



Containing latex: Not manufactured with natural rubber latex

Manufacturer

ACE Medical Devices PVT. LTD.

Distributor & Importer

COLOPLAST A/S

CE Class according to MDD (EU) 93/42/EEC as amended by 2007/47/EC

Is

Indications

Intended to be used for evacuating irrigation solutions, blood clots, tissue collection and/or straining functions during prostate surgery or bladder surgery.

Instructions for use I-BE-COM-CP-00

Sterilisation

Sterile and single use

Ethylene oxide (EO)

Shelf life

3 years

Packaging

Use unit	1 bladder evacuator device in a peel pouch
Order unit	10 use units in a box: 360 x 270 x 190 mm

Change log

Version number	Initials Issue date	Description
1.0	ITPG March 2023	New document creation