

EasiVac

— Evacuator til fjernelse af vævsrester efter prostata resektion



E80x++ EASIVAC® Bladder evacuation device for single use

Manufacturer

COLOPLAST A/S

Distributor

COLOPLAST A/S

CE Class

Is

SterilisationSterile and single use
Ethylene oxyde (EO)**Indications**

Removal of prostate chips after TURP

*Instructions for use SH2096***Shelf life**

5 years

References

Connector for cystoscope	References
WOLF	E80103
STORZ	E80203
OLYMPUS	E80303
ACMI	E80403

DescriptifReservoir capacity (ml) : 170
Filter and anti-reflux valve**Materials** : Connection tube in silicone, reservoir, cap and filter in polypropylène, pump in P.V.C, connector in ABS**Containing latex** : No**Packaging**

Use unit	1 evacuation device in a double peel pouch
Order unit	20 use units in a carton box

Marketing

Regulatory Affairs



Francis Guitteny
International Brand Manager
Porgès - Coloplast Division

4/01/2017

Stéphane Bouché
Regulatory Affairs Manager
Porgès - Coloplast

11/01/2017

