

EasiVac

12. February 2018
Version 3.0

Safety Data Sheet

Based on template version 4.0

Identification of the substance/mixture and of the company/undertaking

Product name: EasiVac
Product code: E80103, E80203, E80303, E80403, E80503
Product information: Disposable Surgical Urology (Urology & Continence)
Manufacturer: Coloplast A/S
Holtedam 1
DK-3050 Humlebaek
Denmark
Telephone +45 49111111
msds@coloplast.com

USA

Coloplast Corp.
1601 West River Road N
Minneapolis, MN 55411
Telephone: +1-800-533-0464
www.us.coloplast.com

Canada

Coloplast Canada Corporation
3300 Ridgeway Drive, Unit 12
Mississauga, Ont. L5L 5Z9
Telephone: +1-877-820-7008
www.coloplast.ca

Europe

Coloplast A/S
Holtedam 1
DK-3050 Humlebaek
Telephone: +45 49 11 11 11
www.coloplast.com

Hazards identification

This product consists primarily of polymer materials. The products pose no immediate hazard.

Composition/information on ingredients

This product does not contain substances classified as hazardous under EC Regulation No. 1272/2008/EC, Annex VI (EU) or OSHA regulation Hazard communication Standard 2012 (USA) in concentrations above 0.1 % (w/w). Main ingredients and packaging materials are listed below.

Chemical name	CAS no.
Polypropylene	9010-79-1
Polyethylene	-
Silicone	-
ABS	100-41-4 107-13-1
PVC (polyvinyl chloride)	-

EasiVac

Packaging:

PA (polyamide)	-
PE (polyethylene)	-
Paper	-
Cardboard	-

Disposal considerations

Dispose the device according to the recommended disposal technology at any approved facility. Under normal private use the product may be disposed of together with other household waste unless local regulations set special requirements. The disposal should always be in compliance with national, federal, state and local regulations. The product should not be discharged to the environment.

US

This product does not meet the criteria for hazardous waste as defined under the Resource Conservation and Recovery Act (RCRA) 40 CFR 261. Under normal private use the product may be disposed of together with other household waste per RCRA 40 RFT 261.4.B1.

European Union

Per The European Waste Catalogue (EWC), in accordance with EC Directive 75/422/ECC, the following Waste Code can be used: 18 01 04 00 wastes whose collection and disposal is not subject to special requirements in view of the prevention of infection (e.g. dressings, plaster casts, linen, disposable clothing, diapers). However, if the waste in view of the prevention of infection needs special requirements, other Waste Codes should be used. Under normal private use the product may be disposed of together with other household waste unless local regulations set special requirements.

Handling and storage

Handling:	See instruction for use
Storage:	Store until use as supplied and at room temperature unless other information is stated on the packaging or on the leaflet.

Other information

This MSDS is supplied as an additional service to the customer. The product is a medical device, which is regulated under the Council Directive 93/42/ECC, Medical Device Directive. The product has been evaluated according to the requirements of medical devices. According to current knowledge this product is considered non-toxic. For further information please contact Coloplast A/S.

Signature Page for VV-0097216 v2.0

Approved by	PLMITI Milena Titoni Junior Biosafety Specialist Technical / Specialist 12-Feb-2018 13:02:43 GMT+0000
-------------	--

Approved by	DKPLYM Peter Lyng Madsen Senior Biosafety Specialist Technical / Specialist 13-Mar-2018 14:20:55 GMT+0000
-------------	--

Signature Page for VV-0097216 v2.0

Document Owner: DKPLYM Peter Lyng Madsen