


# Ureter dilator – dilation af ureter



## Ureteral dilators

### Manufacturer

COLOPLAST A/S

### Distributor

COLOPLAST A/S

### CE class

I (Rule 5)

### Sterilisation

Sterile (Ethylene oxide) and single use

### Indications

Dilatation of the ureter during ureteroscopic procedures.

### Contraindications

Same as for ureteral dilatation.  
Any contraindication to endoscopy.  
Any known allergies to the medical device materials.

*Instructions for use SH2193*

Marketing



### References

CH/FR	Order reference
8-10	RBD010
12-14	RBD014

### Description

Length (cm): 48  
Compatible with a 0.038 inch guidewire  
Distal tip: tapered open

**Material:** Poly Ether Bloc Amide (PEBA)

**Containing latex:** No

### Shelf life

5 years

### Packaging

**Use unit** 1 component per peel pouch

**Order unit** 10 use units in a carton box (mm)  
40 x 100 x 620

Regulatory Affairs



Natacha MIANO  
Regulatory Affairs Manager  
Porgès - a Coloplast division